



1803C Healing Waters

1803C- Health Check Single Top Examination Couch

- 01 Overall approximate dimension: 1900 mm x 560 mm x 840 mm, 75" x 22" x 33" (L x W x H)*
- 02 Overall Width with BP tray: 945 mm, 37"*
- 03 Overall Width with step stool: 970 mm, 38"*
- 04 Mattress platform size: 1900 mm x 565 mm x 65 mm, 75" x 22" x 2.5" (L x W x H)*
- 05 Upper section of box size: 1220 mm x 460 mm x 630 mm, 48" x 18" x 25" (L x W x H)*
- 06 Two side cabinet size: 430 mm x 440 mm x 430 mm, 17" x 17" x 17" (L x W x H)*
- 07 One center cabinet size: 350 mm x 430 mm x 370 mm, 14" x 17" x 15" (L x W x H)*
- 08 BP tray size:
350 mm x 125 mm x 20 mm, 14" x 5" x 1" (L x W x H)*
- 09 Fixed upholstered top in two pieces of 65 mm, 2.5"(W)*
- 10 Headrest adjustable on gas spring
- 11 Upper section of box with three drawers
- 12 Lower section comprises of three cabinets with separate doors and lock on two side doors
- 13 BP apparatus tray provided near headrest
- 14 Inbuilt sliding step stool
- 15 Pre-treated and epoxy powder coated
- 16 Supplied in Semi Knocked Down condition (SKD)
- 17 Safe working load of 170 kgs and patient load bearing capacity of 135 kgs
- 18 Upholstery Colour: Healing Waters and Robust Brown



1803C Robust Brown

Optional Accessories

- 01 Paper roll holder (0175)

1803L



1803L Healing Waters

1803L- Health Check Single Top Examination Couch

- 01 Overall approximate dimension: 1830 mm x 560 mm x 840 mm, 72" x 22" x 33" (L x W x H)*
- 02 Mattress platform size: 1830 mm x 560 mm x 65 mm, 75" x 22" x 2.5" (L x W x H)*
- 03 Upper section of box size: 1220 mm x 460 mm x 630 mm, 48" x 18" x 25" (L x W x H)*
- 04 Two side cabinet size: 340 mm x 430 mm x 430 mm, 13" x 17" x 17" (L x W x H)*
- 05 One center cabinet size: 350 mm x 430 mm x 370 mm, 14" x 17" x 15" (L x W x H)*
- 06 BP tray size: 350 mm x 125 mm x 20 mm, 14" x 5" x 1" (L x W x H)*
- 07 Single section fixed upholstered top
- 08 MS Powder coated drawers and cabinet doors and sliding inbuilt step stool
- 09 Upper section of box with three drawers
- 10 Lower sections comprises of three cabinets with separate doors and lock on two side doors
- 11 BP apparatus tray provided near headrest
- 12 Pre-treated and epoxy powder coated
- 13 Supplied in Semi Knocked Down condition (SKD)
- 14 Upholstery Colour: Healing Waters and Robust Brown



1803L Robust Brown

Optional Accessories

- 01 Paper roll holder (0175)



Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com