



1803M Healing waters

1803M- Health Check Single Top Examination Couch With "U" CUT

- 01 Overall approximate dimension: 1900 mm x 570 mm x 840 mm, 75" x 22" x 33" (L x W x H)*
- 02 Overall Width with BP tray: 920 mm, 36"*
- 03 Mattress platform size: 1840 mm x 570 mm x 65 mm, 72" x 22" x 3" (L x W x H)*
- 04 Upper section of box size: 1220 mm x 460 mm x 630 mm, 48" x 18" x 25" (L x W x H)*
- 05 Two side cabinet size: 430 mm x 320 mm x 430 mm, 17" x 13" x 17" (W x H x D)*
- 06 One center cabinet size: 430 mm x 350 mm x 370 mm, 17" x 14" x 15" (W x H x D)*
- 07 BP apparatus tray size: 356 mm x 127 mm x 20 mm, 14" x 5" x 1" (L x W x H)*
- 08 Single section upholstered top, with "U" cut
- 09 A pair of synthetic rubber padded SS lithotomy poles with rexine ankle straps
- 10 MS powder coated drawers and cabinet doors and sliding inbuilt step stool
- 11 Upper section of box with three drawers
- 12 Lower section comprises of three cabinets with separate doors and lock on two side doors.
- 13 Pre-treated and epoxy powder coated
- 14 Supplied in SKD condition
- 15 Safe working load of 170 kgs and patient weight bearing capacity of 135 kgs
- 16 Upholstery Colour: Healing Waters and Robust Brown



1803M Robust brown

Optional Accessories

- 01 Paper roll holder (0175)

1808



1803M Healing waters

1808- Health Check Examination-Cum-Gynaec Table with Upholstered Top & "U" Cut

- 01 Overall approximate dimension:
1890 mm x 610 mm x 825 mm, 74" x 24" x 33" (L x W x H)*
- 02 Two sections with fixed upholstered tops
- 03 Headrest adjusted on ratchet
- 04 Perineal recess and SS bowl at leg end
- 05 Pair of synthetic rubber padded stainless steel lithotomy poles with rexine ankle straps
- 06 Legs fixed with rubber shoes
- 07 Pre-treated and epoxy powder coated
- 08 Supplied in KDC
- 09 Safe working load of 140 Kgs and patient weight bearing capacity of 135 Kgs
- 10 Upholstery Colour: Healing Waters and Robust Brown



1808 Robust Brown

 *± 10 mm Engineering Variation

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

1803M 1808:05:05.11.23