



1807 Healing Waters

1807- Health Check Examination Table with Upholstered Top

- 01 Overall approximate dimension: 1835 mm x 615 mm x 825 mm, 72"x 24"x 32" (L x W x H)*
- 02 Mild steel tubular construction, machine pressed double bend mild steel sheet
- 03 Two sections with fixed upholstery
- 04 Headrest adjusted on ratchet with four elevation angles
- 05 Legs fixed with rubber shoes
- 06 Pre-treated and epoxy powder coated
- 07 Supplied in KDC
- 08 Safe working load of 140 Kgs and patient weight bearing capacity of 135 Kgs
- 09 Upholstery Colour: Healing Waters and Robust Brown



1807 Robust Brown

*± 10 mm Engineering Variation



*The product images shown are for illustrative purpose only and not be an exact representation of the product. Accessories shown in the product image may not be part of the standard configuration.

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

1807:07:07.11.23

We are interested in growing relationship

