

Motorized xamination Couch



JX2000

Model 1 with
Handset Control



Model 2 with
Foot Switch Control



*The product images shown are for illustrative purpose only and not be an exact representation of the product. Accessories shown in the product image may not be part of the standard configuration.

We are interested in growing relationship





The above image is for representation purpose only. The actual product is available either with foot switch or handset. Safety Side railing (single piece), Adjustable BP Apparatus Tray are part of optional accessories.

Technical Specifications

Platform Dimension*	2010 mm x 655 mm, 79"x26" (LxW)
Minimum Height (Electrically Operated)*	500 mm, 19" without Upholstery + 65 mm Upholstery, 3"
Maximum Height (Electrically Operated)*	860 mm, 34" without Upholstery + 65 mm Upholstery, 3"
Backrest**	0° - 60°
Safe Working Load	170 Kgs
Patient Weight Bearing Capacity^	150 Kgs

Standard Features

- Gas spring assisted backrest
- Handset control or foot switch control
- Battery backup
- Integrated paper roll holder
- User friendly dual side operated utility storage space

Electrical Parameters

Power In (Operating Range)	85 VAC to 265 VAC, 50/60 Hz, Max 5A
Liquid Ingress Protection	IPX6
Electrical Shock Protection	Class 1
Degree of Shock Protection	Type B
Duty Cycle	10 % (Two minutes for every eighteen minutes)

Safety Standards

General safety for medical electrical equipment	IEC 60601-1-1
EMI/EMC	IEC 60601-1-2
General safety for medical equipment	IEC 60601-2-52

Variants

Features	Model 1	Model 2
Handset	✓	✓
Foot Switch Control	✓	✓

*± 10 mm Engineering Variation ** ± 5° Engineering Variation ^ When couch is configured with standard accessories

Add-On Accessories



MS Safety Side Railing
(Single piece can be installed on either side of the couch)



MS Adjustable BP Apparatus Tray

Upholstery Colour Options



Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

JX2000:09:22.11.23



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark 'Midmark India' and the logo. Reproduction of this catalogue in any form is liable for legal action.