



## Motorized Gynaecology Examination Couch



JX4000

### Clinical Positions



## Technical Specifications

Overall Length* (Inclusive of foldable leg section of 490 mm)	1695 mm, 67"
Overall Approximate Width*	920 mm, 36"
Upholstered Platform Width*	680 mm, 27"
Upholstery Thickness*	65 mm, 3"
Maximum Height (Electrically Operated)*	920 mm, 36"
Minimum Height (Electrically Operated)*	560 mm, 22"
Backrest Angle**	0° - 70°
Safe Working Load	180 Kgs
Patient Weight Bearing Capacity^	175 Kgs

## Standard Features

- Motorised backrest and height adjustment
- Wired handset control for patient and attendant for backrest and height adjustment
- Manually operated pelvic tilt
- Pair of Lithotomy crutches with cover
- Integrated paper roll holder
- Removable pan
- Handgrips
- Battery backup


## Electrical Specifications

Power In (Operating Range)	85 VAC to 265 VAC, 50/60 Hz, Max 5A
Liquid Ingress Protection	IPX6
Electrical Shock Protection	Class 1
Degree of Shock Protection	Type B
Duty Cycle	10 % (two minutes for every eighteen minutes)

## Safety Standards

General safety for medical electrical equipment	IEC 60601-1-1
EMI/EMC	IEC 60601-1-2
General safety for medical equipment	IEC 60601-2-52

## Overall Dimension



## Upholstery Colour Options



Moss



Navy



Raisin Grape



Shadow



Pista



Healing Waters



Robust Brown

### Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070  
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

JX4000:06:05.10.23



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.