



9004A

New Modular Crash Cart 9004A

Designed to attend patients instantly with life saving equipment and drugs.
The Crash Cart 9004A is super-efficient and user-friendly & convenient during emergencies.

We are interested in growing relationship



Overall Approximate Dimension

Pushing Handle to Oxygen Cylinder Cage Provision *	1020 mm x 680 mm, 40" x 27" (L x W)
Platform Dimension: (Top Polymer Moulded Panel) *	560 mm x 370 mm, 22" x 15" (L x W)
Buffer to Buffer Dimension *	820 mm x 680 mm, 32" x 27" (L x W)
Height (Floor to Top Polymer Moulded Panel) *	950 mm, 37" (H)
Maximum Height (Floor to Eye-level Detachable Fully SS Top Shelf) *	1645 mm, 65" (H)

Technical Specifications

Cardiac Massage Board (Dimensions) *	710 mm x 390 mm, 28" x 15" (L x W) 8 mm Thickness
First Drawer ABS Tray *	578 mm x 360 mm x 65 mm, 23" x 15" x 2.5" (L x W x H)
Second & Third Drawer ABS Tray *	578 mm x 360 mm x 145 mm, 23" x 15" x 5.71" (L x W x H)
Forth Drawer ABS Tray *	578 mm x 360 mm x 225 mm, 23" x 15" x 8.85" (L x W x H)
Provision for De-Fib Tray *	335 mm x 235 mm x 20 mm, 13" x 9" x 0.8" (L x W x H)

Standard Features

<ul style="list-style-type: none"> Four drawers with central locking
<ul style="list-style-type: none"> Top & second drawer each will have poly-carbonate partition plates of 1 set includes (3 long plates + 3 short plates)
<ul style="list-style-type: none"> File / chart holder
<ul style="list-style-type: none"> Trash can (1 no.)
<ul style="list-style-type: none"> Catheter holder (2 nos.)
<ul style="list-style-type: none"> Oxygen cylinder cage
<ul style="list-style-type: none"> Two shelves, one ss & one with multicoloured removable bins at eye-level along with IV pole with single hook
<ul style="list-style-type: none"> Twin wheel high grade synthetic castors 125 mm (5") dia castors two with brake and two without brake buffers at 4 corners

Add-On Accessories

<ul style="list-style-type: none"> De-fib tray holder platform of 310 mm x 210 mm, 12" x 8" (L x W) (0130H) *
<ul style="list-style-type: none"> Two shelves, one ss & one with multicoloured removable bins at eye-level along with iv rod with single hook" (0130G)
<ul style="list-style-type: none"> Scissor holder (0130C)
<ul style="list-style-type: none"> Additional set of polycarbonate partition for third drawer (0130D)
<ul style="list-style-type: none"> Additional set of polycarbonate partition plates for 3rd middle single drawer (0130E)
<ul style="list-style-type: none"> Set of polycarbonate partition plates for fourth drawer (0130F)

* Engineering Variation ± 10 mm

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

9004A-06:03.12.23



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.