


# Emergency Drug Cart



## ICU Solutions Emergency Drug Cart

Designed to attend patients instantly with life saving equipment and drugs. The emergency drug cart is super-efficient and user-friendly & convenient during emergencies.

Backside View



We are interested in growing relationship



## 9005 - Technical Specifications

Overall Approximate Dimension	
Dustbin to Opposite Side Basket*	860 mm x 530 mm (LxW)
Platform Dimension* (Top Polymer Moulded Panel) (LxW)	550 mm x 400 mm (LxW)
Height* (Floor to Top Polymer Moulded Panel)	980 mm (H)
Sliding Slide Tray* (LxW)	355 mm x 320 mm (LxW)
IV Pole Height and Extension* (From Top Platform)	225 mm to 760 mm
Accessories	
Cardiac Massage Board* (LxW)	525 mm x 390 mm (LxW)
Defib Tray* (LxWxH)	407 mm x 372 mm x 15 mm (LxWxH)
Dustbin (2 nos.)	Yes
Provision for Oxygen Cylinder	Yes
Five Drawers with Central Locking	Yes (Interior with easily adjustable section)
Basket	Yes
Tray Basket	Yes
Medical Waste Bin with Central Locking	Yes
Provision for Defibrillator	Yes
Label Cards for Drawers	Yes
High Grade Castors (Two with brake and two without brake)	Yes

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.

\* Engineering Variation ± 10 mm



### Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070  
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

