


CX4000 Plus



CX4000 Plus

Presenting the all new CX4000 Plus, a thoughtfully crafted variant designed for the challenging healthcare environments, keeping caregiver's needs in mind.

It is packed with the full set of clinical positions like height, backrest, trendelenburg, reverse trendelenburg, knee rest on individual smooth crank systems and a multi position DVT end platform, making it a perfect choice for critical care and high dependency areas. Armed with an elegantly designed and smoothly operated and split type polymer molded safety side rails, it not only offers safety and hygiene but also brings treatment areas to life. The bed is engineered for trouble free service life with a tough frame design, strong support linkages and high grade synthetic castors.



*The product images shown are for illustrative purpose only and not be an exact representation of the product. Accessories shown in the product image may not be part of the standard configuration.

We are interested in growing relationship



Technical Specifications

Buffer to Buffer Dimension*	2140 mm X 1020 mm, 85" X 40" (L X W)
Platform Dimension*	2080 mm X 955 mm, 82" X 38" (L X W)
Mattress Platform*	1890 mm X 860 mm, 76" X 34" (L X W)
Minimum Height * (without mattress)	530 mm, 21"
Maximum Height * (without mattress)	770 mm, 30"
Backrest**	0° - 70°
Knee Rest**	0° - 45°
Trendelenburg**	0° - 8°
Reverse Trendelenburg**	0° - 8°
Safe Working Load	150 kgs
Patient Load Bearing Capacity^	110 kgs

Standard Features

• Raised backrest
• Raised upper leg
• Height adjustment
• Trendelenburg/ reverse trendelenburg
• DVT end platform/ calf elevation
• Four section metallic perforated top
• Pre-treated and epoxy powder coated MS parts
• 125 mm high grade synthetic castors, two with brake and two without brake
• Polymer moulded split type safety side rails
• Polymer moulded head and foot boards (pair)
• Safety buffers at four corners
• Removable handle
• Urine bag holder

Variants

Features	CX4000 Plus 380	CX4000 Plus 360
Coverage span	Full length	¾ th length ⁵
Angle indicator for backrest, upper leg, trendelenburg/ reverse trendelenburg	Available	-

Add-On Accessories

• Single section premium mattress	• SS heavy duty IV pole with four hooks
• Utility IV pole with two hooks	• Urine Bag Holder
• Utility IV pole with four hooks	• Four section mattress
• Telescopic IV pole with two hooks	• Single section mattress
• Telescopic IV pole with four hooks	• Vertical oxygen cylinder cage




Oxygen Cylinder Cage




Orthopaedic Traction Pulley



Stainless Steel Lifting Pole



Moulded Chart Holder



3rd Set of Safety Side Railings

- All variants are available with X-Ray permeable backrest (XR) feature as an add-on accessory, which are factory fitted.

* ±10 mm Engineering variation. ** ± 5° Engineering variation. ^When bed is configured with standard accessories.

⁵ Can be converted into full length by adding 3rd set railing.

Note: Option of fixed polymer moulded handle available in lieu of removable handle.

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

CX4000 Plus: 09:13.03.24

