

ICU Solutions



JANAK[®]
Because we care.



CX4000A Comfort Bed

This bed combines the full features of a hospital bed that comprise of back adjustment, upper leg, height and trendelenburg / reverse trendelenburg on four separate crank systems. This bed is beneficial to clinics, nursing homes, and hospitals. This adjustable manual hospital bed is designed to deliver Midmark India quality and safety at an affordable price, providing comfort to your patients and caregivers.

*The product images shown are for illustrative purpose only and not be an exact representation of the product. Accessories shown in the product image may not be part of the standard configuration.



ICU Solutions

Technical Specifications

Buffer to Buffer Dimension*	2140 mm x 970 mm, 84" x 38" (L x W)
Platform Dimension*	2090 mm x 920 mm, 82" x 36" (L x W)
Mattress Platform*	1950 mm x 820 mm, 77" x 32" (L x W)
Buffer to Handle*	2325 mm, 92"
Minimum Height* (without mattress)	500 mm, 20"
Maximum Height* (without mattress)	740 mm, 29"
Backrest**	0° - 70°
Knee Rest**	0° - 50°
Trendelenburg**	0° - 9°
Reverse Trendelenburg**	0° - 9°
Safe Working Load	150 kgs
Patient Load Bearing Capacity^	110 kgs

Add-On Accessories

- Four section mattress with 4" thick PU foam
- Telescopic IV pole with 2 hooks
- Telescopic IV pole with 4 hooks
- SS heavy duty IV pole with 2 hooks
- SS heavy duty IV pole with 4 hooks
- Oxygen cylinder cage (epoxy coated)
- Orthopedic traction pulley
- SS lifting pole
- Fixed polymer moulded handle (set of four)
- Removable handle
- Moulded chart holder
- Urine bag holder

* ±10 mm Engineering variation. ** ± 5° Engineering variation. ^When bed is configured with standard accessories
Note: Option of fixed polymer moulded handle available in lieu of removable handle.

Standard Features

- Raised backrest
- Raised kneerest
- Trendelenburg / reverse Trendelenburg
- DVT end platform / calf elevation
- Four section metallic perforated top
- Pre-treated and epoxy powder coated MS parts
- 125 mm high grade synthetic castors, two with brake and two without brake
- Buffers at 4 corners for safety
- Provision for IV Pole
- Removable handle
- Easily removable head and foot boards

Safety Standard

General safety for medical beds IEC 60601-2-52

Railings

- Collapsible type full-length type SS safety side railing with epoxy powder coated frame (pair)
- Swing down type full-length epoxy powder coated MS safety side railing (Pair) (factory fitted Only)
- Collapsible type epoxy powder coated vertical safety side railing (pair)
- Collapsible type epoxy powder coated safety side railing (pair)
- Collapsible type epoxy powder coated MS safety side railing with SS tubes (pair)

Head and Foot Boards

- Polymer moulded head and foot boards (pair)
- SS laminated head & foot boards (pair)



Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

We are interested in growing relationships



CX4000A: 07:19.05.23