

# ICU Solutions



**JANAK**<sup>®</sup>  
Because we care.



## CX4000B Comfort Bed

The comfort bed consists of fixed height, four sections with mild steel perforated sheet top. It comprises of backrest, trendelenburg, and reverse trendelenburg position on the crank mechanism with a removable handle. Easily removable head and footboards and high-grade synthetic body castors of 125 mm diameter, two with brake and two without brake.

\*The product images shown are for illustrative purpose only and not be an exact representation of the product. Accessories shown in the product image may not be part of the standard configuration.



# ICU Solutions

## Technical Specifications

Buffer to Buffer Dimension*	2140 mm x 970 mm, 84" x 38" (L x W)
Platform Dimension*	2090 mm x 920 mm, 82" x 36" (L x W)
Mattress Platform*	1960 mm x 820 mm, 77" x 32" (L x W)
Buffer to Handle*	2260 mm, 89"
Fixed Height*(without mattress)	600 mm, 24"
Backrest**	0° - 70°
Knee Rest**	0° - 50°
Trendelenburg**	0° - 10°
Reverse Trendelenburg**	0° - 7°
Safe Working Load	150 kgs
Patient Load Bearing Capacity^	110 kgs

## Standard Features

- Raised backrest
- Raised kneerest
- Trendelenburg / reverse Trendelenburg
- DVT end platform/ calf elevation
- Four section metallic perforated top
- Pre-treated and epoxy powder coated MS parts
- 125 mm high grade synthetic castors, two with brake and two without brake
- Buffers at 4 corners for safety
- Provision for IV Pole
- Removable handle
- Easily removable head and foot boards

## Safety Standard

General safety for medical beds

IEC 60601-2-52

## Add-On Accessories

- Four section mattress with 4" thick PU foam
- Telescopic IV pole with 2 hooks
- Telescopic IV pole with 4 hooks
- Oxygen cylinder cage (epoxy coated)
- Orthopedic traction pulley
- SS lifting pole
- Fixed polymer moulded handle (set of three)
- Removable handle
- Moulded chart holder
- Urine bag holder

## Railings

- Collapsible type full-length type SS safety side railing with epoxy powder coated frame (pair)
- Swing down type full-length epoxy powder coated MS safety side railing (Pair) (factory fitted Only)
- Collapsible type epoxy powder coated vertical safety side railing (pair)
- Collapsible type epoxy powder coated safety side railing (pair)
- Collapsible type epoxy powder coated MS safety side railing with SS tubes (pair)

## Head and Foot Boards

- Polymer moulded head and foot boards (pair)
- SS laminated head & foot boards (pair)

\* ±10 mm Engineering variation. \*\* ±5° Engineering variation. ^When bed is configured with standard accessories  
Note: Option of fixed polymer moulded handle available in lieu of removable handle.



## Midmark (India) Private Limited

**Head Office:** Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070  
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

We are interested in growing relationships



CX40008: 04:19:05.23