



ICU Solutions

DX4305 4DX205



Our DX series are the beds for all requirements. Modular construction and wide range of features make DX the all-round bed.

From standard to optional features and a blend of efficiency and durability that's how DX makes an impression.



*The product images shown are for illustrative purpose only and not be an exact representation of the product. Accessories shown in the product image may not be part of the standard configuration.

We are interested in growing relationship



Technical Specifications

Overall Dimension (Buffer to Buffer)*	2150 mm L x 1025 mm W
Platform Frame Dimension*	2080 mm L x 955 mm W
Buffer to Handle*	2315 mm
Mattress Platform*	1910 mm L to 865 mm W
Mattress Platform Height*	530 mm to 780 mm
Backrest**	0°-70°
Safe Working Load	200 Kgs
Patient Load bearing capacity^	175 Kgs

Electrical Parameters

Power In (Mains Power Supply)	100 - 240 V ~, 50/60 Hz., Max 5A
Electrical Shock Protection	Class I
Degree of Shock Protection	Type B
Liquid Ingress Protection	IPX6
Duty Cycle	10 % (Two minutes for every eighteen minutes)

Safety Standards

General safety for medical beds	IEC 60601-2-52
General safety for medical electrical equipment	IEC 60601-1-1
EMI/EMC	IEC 60601-1-2

Variants

Features	DX4305 0392	DX4305 380	DX4305 360	DX4205 0392	DX4205 380	DX4205 360
Knee Rest**	0°- 45°	0°- 42°	0°- 42°	0°- 45°	0°- 45°	0°- 45°
Upper Leg	Motorized	Motorized	Motorized	Manual	Manual	Manual
Coverage span	Full length	Full length	¾ th length ⁵	Full length	Full length	¾ th length ⁵
Screw mechanism for trendelenburg & rev. trendelenburg tilt adjustment with single removable handle to operate the above mentioned contour position	Available	Available	Available	-	-	-
Screw mechanism for upper leg, trendelenburg & rev. trendelenburg tilt adjustment with single removable handle to operate the above mentioned contour position	-	-	-	Available	Available	Available
Trendelenburg	0°- 10°	0°- 10°	0°- 10°	0°- 10°	0°- 9°	0°- 9°
Reverse Trendelenburg	0°- 7°	0°- 7°	0°- 7°	0°- 7°	0°- 9°	0°- 9°

Add-On Accessories



Oxygen Cylinder Cage



Orthopaedic Traction Pulley



Stainless Steel Lifting Pole



Moulded Chart Holder



3rd Set of Safety Side Railings

- Provision for Heavy Duty and Telescopic IV Pole
- Single Section Premium Mattress and Four Section PU Foam Mattress
- Battery backup
- All variants are available with X-Ray permeable backrest (XR) feature as an add-on accessory, which are factory fitted.

*±10 mm Engineering variation. **± 5° Engineering variation. ^When bed is configured with standard accessories.

⁵Can be converted into full length by adding 3rd set railing.

Note: Option of fixed polymer moulded handle available in lieu of removable handle.

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

DX4305 4DX205: 11:24.01.24

Standard Features

- Motorized backrest & height adjustment
- Wired handset control
- Dual side manually operated CPR lever
- Calf elevation/DVT platform
- Angle indicator for backrest, trendelenburg /reverse trendelenburg and Quick Assist function for DX4305
- Polymer moulded split type safety side rails
- Polymer moulded head & foot panels, easily removable without using any tools
- Four corner safety rubber buffers
- 125 mm. diameter single wheel high grade synthetic body castors with accessible foot brake mechanism on two castors
- Inbuilt SMPS
- Easily removable polypropylene 5 section mattress top
- Single side urine bag holder
- Provision for SS heavy duty IV pole at 4 corners
- All MS parts are pre-treated and powder coated
- Removable Handle

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.