



stanchion

PX4005

Stanchion, a modern hospital bed, is uniquely designed to provide further comfort to the patients and make the nursing staff work easier. Stanchion offers quality beyond comparison blended with a wide range of features such as patient comfort, caregiver convenience, ease of mobility, column positions ensure excellent stability and is engineered to international standards. This new bed from Midmark India is focused on helping physicians, and care teams deliver the most effective patient care by keeping its commitment to safety and quality.



\*The product images shown are for illustrative purpose only and not be an exact representation of the product. Accessories shown in the above image are not in standard configuration as 3rd set is not present in the default configuration.

We are interested in growing relationship



## Technical Specifications

Overall Dimension*	2150 mm L x 1020 mm W
Mattress Platform*	1920 mm L x 860 mm W
Mattress Platform Height*	460 mm to 790 mm (Without Mattress)
Bed Extension*	215 mm
Backrest**	0°-70°
Knee Rest**	0°-45°
Trendelenburg**	0°-13°
Reverse Trendelenburg**	0°-13°
Safe Working Load	165 Kgs
Patient Load Bearing Capacity^	200 Kgs

## Electrical Parameters

Power In (Mains Power Supply)	100 - 240 V ~, 50/60 Hz., Max 5A
Electrical Shock Protection	Class I
Degree of Shock Protection	Type B
Liquid Ingress Protection	IPX6
Duty Cycle	10 % (Two minutes for every eighteen minutes)

## Safety Standards

General safety for medical beds	IEC 60601-2-52
General safety for medical electrical equipment	IEC 60601-1-1
EMI/EMC	IEC 60601-1-2

## X-Ray Variant

X-Ray Permeable Backrest & Cassette Holder


- Stanchion XR

## Standard Features



Embedded panel with lock button

Attendant control panel with lock button



Column Based Architecture

- One touch key for cardiac chair, auto CPR, emergency head down
- Motorized backrest, knee rest, height adjustment, auto contour, trendelenburg & reverse trendelenburg
- Dual side manually operated CPR lever
- Calf elevation/DVT platform
- Polymer moulded split type safety side rails (2 sets)
- Easily removable polymer moulded head and foot panels.
- Four corner safety rubber buffers of 100 mm diameter
- 125 mm diameter twin wheel high-grade synthetic body castors with central locking mechanism
- Inbuilt SMPS
- Integrated urine bag holder
- Built-in bed extension
- Single section premium mattress
- Mattress arrester (6 nos.)
- Under bed light

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.

\* ±10 mm Engineering variation \*\* ± 5° Engineering variation ^When bed is configured with standard accessories. #Factory fitted.

## Add-On Accessories



Oxygen Cylinder Cage



Orthopaedic Traction Pulley



Stainless Steel Lifting Pole



Moulded Chart Holder



Bed Extension#

- Four Section 4" PU Foam Mattress
- Heavy Duty IV Pole with Two Hooks / Four Hooks
- Bed Extension Mattress
- Battery backup
- Urine Bag Holder
- Vertical Oxygen Cylinder Cage

### Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070  
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

Stanchion: 07:07.01.24

