



0104 – Arrive – Paediatric Bed with Drop Side Rails

- Overall approximate dimension: 1490 mm x 830 mm x 600 mm, 59" x 33" x 24" (L x W x H)*
- Bed frame approximate dimension: 1370 mm x 760 mm, 54" x 30" (L x W)*
- Perforated MS sheet top
- Both bows of 1065 mm (42") H, with one horizontal, and seven vertical supports and rubber feet*
- Full length drop side rails
- Provision for mosquito net poles
- Pre-treated and powder coated
- Supplied in KDC

Optional Accessories

Rexine covered mattress of 100 mm thickness (0111D)

*± 10mm engineering variation

CE

Maternal & Gynaec Solutions



1083 – Crib Carrier

- Overall approximate dimension: 920 mm x 460 mm x 950 mm, 36" x 18" x 37" (L x W x H)*
- Approximate crib size: 760 mm x 375 mm x 305 mm, 30" x 15" x 12" (L x W x H)*
- 3" PU Foam Mattress with PVC Rexine
- Tubular frame fitted with swiveling castors of 125 mm diameter synthetic body castors two with brake and two without brake
- Pre-treated and epoxy powder coated


*± 10mm engineering variation



Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

We are interested in growing relationships



0104-1083:04:20.04.21