



1504A

1504A - Obstetric Labour Table (MS & with SS Top)

- 01** Overall approximate extended dimension: 1880 mm x 910 mm x 800 mm, 74" x 36" x 31" (L x W x H)*
- 02** Mattress platform dimension: 1805 mm x 655 mm, 71" x 26" (L x W)*
- 03** Backrest adjustable on ratchet: 745 mm x 625 mm, 29" x 25" (L x W)*
- 04** Middle section with "U" cut: 350 mm x 650 mm, 14" x 26" (L x W)*
- 05** Leg section can be telescoped under the main frame completely for lithotomy positions: 700 mm x 600 mm, 28" x 24" (L x W)*
- 06** Size of railing: 685 mm x 195 mm, 27" x 8" (L x H)*
- 07** Backrest**: 0°-75°
- 08** Foot operated gas spring assisted trendelenburg**: 0°-12°
- 09** Complete with stainless steel head & side rails, hand grips, a pair of rubber padded stainless steel lithotomy poles with rexine ankle straps and stainless steel bowl under 'U' cut
- 10** Three section SS top & MS tubular construction
- 11** Supplied in SKD condition

Add on Accessories

- 01** Lithotomy crutches in lieu of lithotomy poles
- 02** Rexine covered mattress of 25 mm, 1" thickness
- 03** Stainless steel IV pole- 2 hooks & 4 hooks



*The product images shown are for illustrative purpose only and not be an exact representation of the product. Accessories shown in the product image may not be part of the standard configuration.

We are interested in growing relationship



1502A



1502A

1502A - Obstetric Labour Table (Fully SS)

- 01 Overall approximate extended dimension: 1880 mm x 910 mm x 800 mm, 74" x 36" x 31" (L x W x H)*
- 02 Mattress platform dimension: 1805 mm x 655 mm, 71" x 26" (L x W)*
- 03 Backrest adjustable on ratchet: 745 mm x 625 mm, 29" x 25" (L x W)*
- 04 Middle section with "U" cut: 350 mm x 655 mm, 14" x 26" (L x W)*
- 05 Leg section can be telescoped under the main frame completely for lithotomy positions: 700 mm x 600 mm, 28" x 24" (L x W)*
- 06 Size of railing: 685 mm x 195 mm, 27" x 8" (L x H)*
- 07 Backrest**: 0°-75°
- 08 Foot operated gas spring assisted trendelenburg**: 0°-12°
- 09 Complete with stainless steel side rails, hand grips, a pair of rubber padded stainless steel lithotomy poles with rexine ankle straps and stainless steel bowl under 'U' cut
- 10 Three section SS top and SS tubular construction
- 11 Supplied in SKD condition

Add on Accessories

- 01 Lithotomy crutches in lieu of lithotomy pole
- 02 Rexine covered mattress of 25 mm, 1" thickness
- 03 Stainless steel IV pole – 2 hooks & 4 hooks

*± 10 mm Engineering Variation and ** ± 5° Engineering Variation

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

1502A-1504A:04:10.02.24