


1526



Lateral Tilt Position



Flex Position







Total Viewing of 'KUB' Region



Trendelenburg & Reverse Trendelenburg

Dimensions



We are interested in growing relationship



Technical Data

Overall Approximate Dimension*	1860 mm x 570 mm (LxW), 73" x 22" without accessories
--------------------------------	---

Electro-Mechanical Adjustments

Height*	730 mm to 930 mm, 29" x 37"
Trendelenburg**	0 ° - 25°
Rev. Trendelenburg**	0 ° - 25 °
Lateral Tilt**	0 ° - 15 °

Manual Adjustments

Backrest Adjustment on Gas Spring Mechanism**	0 °- 65 ° Up/0 °- 40° Down
Detachable Head Rest	Available
Detachable Leg Rest	Available

Electrical Specifications

Power In (Mains Power Supply)	Nominal 230 VAC
	90 VAC to 300 VAC Operating Range with Battery
	45Hz to 55 Hz
Electrical Shock Protection	Class I
Degree of Shock Protection	Type B
Liquid Ingress Protection	IP54
Batterie Backup	Available
Rechargeable Batteries	2 x 12 Volt Sealed Lead Acid
Duty Cycle	10 % (Two minutes for every eighteen minutes)

Standard Accessories

- Four section mattress (1 set)
- Lithotomy crutches with radial clamp (2 nos.)
- Arm rest (2 nos.)
- Anaesthesia screen with clamp (1 no.)
- Body strap (1 no.)
- Perineal cutout filler (1 no.)

Features


- SS Siderails (12 mm x 10 mm) on both sides of the table
- Column covers the electro-mechanical actuator system
- Top made of MS tubular framework and fitted with X-Ray permeable 8 mm thick top
- Table has heavy metal base covered with moulded ABS
- Table is fitted with non rusting 50 mm diameter castors
- Foot lever operated table lock
- All MS parts are pre-treated and epoxy powder coated

Safety Standards

- General safety for medical electrical equipment IEC 60601-1-1
- EMI/EMC IEC 60601-1-2

* ±10 mm Engineering variation. ** ± 2° Engineering variation.

Add-On Accessories




Dual Joint Headrest with Mattress



Shoulder Support with Pad




Lateral Support



Pelvic Extension



TUR Drainage Set



Body Strap

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

OTS-A:06:24.02.24



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.