

Revolving Stool



1601A - Revolving Stool

- Tubular tripod based with SS revolving top with fixed cushion of 65 mm(3") thickness*
- Height adjustable from approximately 550 mm (22") to 760 mm (30") with cushion by accurately machined screw mechanism*
- Legs fitted with PVC shoes
- MS tubular construction
- Pre-treated and powder coated finish



1602A- Revolving Stool Fully Stainless Steel

- Tubular tripod based with SS revolving top with fixed cushion of 65 mm (2") thickness*
- Height adjustable from approximately 550 mm (22") to 760 mm (30") with cushion by accurately machined screw mechanism*
- Legs fitted with PVC shoes
- SS tubular construction



1601 - Revolving Stool

- Tubular tripod based with SS revolving top
- Height adjustable from approximately 460 mm (18") to 665 mm (26") by accurately machine screw mechanism*
- Legs fitted with rubber shoes
- MS tubular construction
- Pre-treated and powder coated finish



1602- Revolving Stool Fully Stainless Steel

- Tubular tripod based with SS revolving top
- Height adjustable from approximately 460 mm (18") to 665 mm (26") by accurately machined screw mechanism*
- Legs fitted with rubber shoes
- SS tubular construction



1811 – Two Tier Bowl Stand

- Overall approximate height: 850 mm, 33"(H)*
- Mild steel tubular frame mounted on three swivel, synthetic body castors of 50 mm (2") diameter, without brake
- Two stainless steel bowls of approximately 375 mm, 15" diameter*
- Pre-treated and powder coated

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.

CE *± 10 mm engineering variation

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 www.midmark.in | india@midmark.com

OTS-E:07:28.02.24

