

Revolving Stool



1605A - Revolving Stool

- Tubular quadrapod based with SS revolving top with fixed cushion of 58 mm (2.5") thickness
- Height adjustable from approximately 520 mm (21") to 730 mm (29") with cushion by accurately machined screw mechanism*
- Legs fitted with PVC shoes
- MS tubular construction
- Pre-treated and powder coated finish



1606A - Revolving Stool Fully Stainless Steel

- Tubular quadrapod based with SS revolving top with fixed cushion of 58 mm(2.5") thickness
- Height adjustable from approximately 520 mm (21") to 730 mm (29") with cushion by accurately machined screw mechanism*
- Legs fitted with PVC shoes
- SS tubular construction

*± 10 mm Engineering Variation



1605 - Revolving Stool

- Tubular quadrapod based with SS revolving top
- Height adjustable from approximately 460 mm (18") to 670 mm (26") with cushion by accurately machined screw mechanism*
- Legs fitted with rubber shoes
- MS tubular construction
- Pre-treated and powder coated finish
- SS ring for footrest



1606 - Revolving Stool Fully Stainless Steel

- Tubular quadrapod based with SS revolving top
- Height adjustable from approximately 460 mm (18") to 670 mm (26") with cushion by accurately machined screw mechanism*
- Legs fitted with rubber shoes
- SS tubular construction
- SS ring for footrest

*±10mm Engineering variation

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 www.midmark.in | india@midmark.com

OTS-A:04:28.02.24



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.