



1812 – Single Tier Bowl Stand

- Overall approximate dimension: 850 mm, 33"(H)*
- Stainless steel bowl size: 375 mm, (15") diameter*
- Twin wheel castors of 50 mm, (2") diameter
- Pre-treated and powder coated



1612 – Stainless Steel Kick Bucket

- Stainless steel tubular framework fitted with three 125 mm, (5") diameter synthetic body castors, without brake
- Stainless steel bowl of approximate 375 mm, (15") diameter*
- Overall height: 330 mm, (15")*
- Supplied in SKD condition



1064 – Stainless Steel Soiled Linen Trolley

- Overall approximate dimension: 960 mm X 510 mm, 38" X 20" (L X W)*
- Stainless steel tubular framework fitted with three 125 mm (5") diameter synthetic body castors, without brake
- Supplied with canvas bag
- Supplied in SKD condition



1633 – Stainless Steel Mayo's Table

- Adjustable height of SS tray from 665 mm (min) to 1160 mm (max), 26" to 46"*
- Tray dimension: 560 mm X 400 mm, 22" X 16" (L X W)*
- SS Mayo's Table with three castors 50mm diameter
- Supplied in SKD condition
- Safe working load: 5 kgs

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.

CE *± 10 mm engineering variation

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
 Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 www.midmark.in | india@midmark.com

1612-1633-1064-1812:07:26.02.24

