



1625



1625A



1625B

1625 - Stainless Steel Instrument Trolley

- Overall approximate dimension: 680 mm x 450 mm x 900 mm, 27"x 18"x 35" (L x W x H)*
- Stainless steel tubular frame mounted on four 125 mm diameter castor wheels with synthetic body, two with brake & two without brake
- Two stainless steel shelves with protective railing on all three sides
- Supplied in Semi Knocked Down condition (SKD)

1625A - Stainless Steel Instrument Trolley

- Overall approximate dimension: 1067 mm x 610 mm x 952 mm, 42"x 24"x 37" (L x W x H)*
- Stainless steel tubular frame mounted on four 125 mm diameter castor wheels with synthetic body, two with brake & two without brake
- Two stainless steel shelves with protective railing on all three sides
- Supplied in Semi Knocked Down condition (SKD)

1625B - Stainless Steel Instrument Trolley

- Overall approximate dimension: 914 mm x 508 mm x 900 mm, 36"x 20"x 35" (L x W x H)*
- Stainless steel tubular frame mounted on four 125 mm diameter castor wheels with synthetic body, two with brake & two without brake.
- Three stainless steel shelves with protective railing on three sides.
- Supplied in Semi Knocked Down condition (SKD)

*± 10 mm engineering variation



OT Solutions



1067



1068



1624

1067 - Stainless Steel Dressing Trolley with SS Bowl & Bucket

- Overall approximate dimension: 1015 mm x 510 mm x 920 mm, 40"x 20"x 36" (L x W x H)*
- Approximate shelf dimension: 755 mm x 500 mm, 30"x 20" (L x W)*
- Stainless steel tubular frame mounted on four 125 mm diameter castor wheels with synthetic body, two with brake & two without brake
- Two stainless steel shelves with protective railing on all four sides
- Supplied in Semi Knocked Down condition (SKD)

1068- Stainless Steel Instrument Trolley

- Overall approximate dimension: 760 mm x 510 mm x 920 mm, 30"x 20"x 36" (L x W x H)*
- Approximate shelf dimension: 755 mm x 500 mm, 30"x 20" (L x W)*
- Stainless steel tubular frame mounted on four 125 mm diameter castor wheels with synthetic body, two with brake & two without brake
- Two stainless steel shelves with protective railing on all four sides
- Supplied in Semi Knocked Down condition (SKD)

1624 - Stainless Steel Instrument Trolley with SS Bowl & Tray


- Out to Out dimension with bowl: 860 mm x 650 mm, 30"x 20" (L x W)*
- Overall approximate dimension: 680 mm x 450 mm x 916 mm, 27"x 18"x 36" (L x W x H)*
- Stainless steel tubular frame mounted on four 125 mm diameter castor wheels with synthetic body, two with brake & two without brake
- Two stainless steel shelves with protective railing on three sides
- With stainless steel bowl & tray
- Supplied in Semi Knocked Down condition (SKD)

*± 10 mm engineering variation



Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com



We are interested in growing relationships

OTS-D:03:20.07.21