



1626 Stainless Steel Instrument Trolley

- Overall approximate dimension: 680 mm x 450 mm x 1000 mm, 27" x 18" x 39" (L x W x H)*
- Stainless steel tubular frame mounted on four 125 mm diameter castor wheels with synthetic body, two with brake & two without brake
- Three stainless steel shelves with protective railing on three sides
- Supplied in Semi Knocked Down condition (SKD)



1626A Stainless Steel Instrument Trolley

- Overall approximate dimension: 1067 mm x 610 mm x 1000 mm, 42" x 24" x 39" (L x W x H)*
- Stainless steel tubular frame mounted on four 125 mm diameter castor wheels with synthetic body, two with brake & two without brake
- Three stainless steel shelves with protective railing on three sides
- Supplied in Semi Knocked Down condition (SKD)



1626B Stainless Steel Instrument Trolley

- Overall approximate dimension: 914 mm x 508 mm x 1000 mm, 36" x 20" x 39" (L x W x H)*
- Stainless steel tubular frame mounted on four 125 mm diameter castor wheels with synthetic body, two with brake & two without brake
- Three stainless steel shelves with protective railing on three sides
- Supplied in Semi Knocked Down condition (SKD)

CE *± 10 mm engineering variation

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

OTS-A:05:26.02.24

We are interested in growing relationship

