

OT Solutions



JANAK[®]
Because we care.



9012



9014

Modular Monitor Trolley



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.



OT Solutions

Technical Specifications

	9012	9014
Overall Approximate Dimension*	680 mm x 670 mm x 1580 mm, 27" x 26" x 62" (L x D x H) *	680 mm x 780 mm x 1360 mm, 27" x 31" x 54" (L x D x H) *
Overall Approximate Dimensions with Gas Cylinder Cage*	680 mm x 670 mm x 1580 mm, 27" x 26" x 62" (L x D x H) *	680 mm x 700 mm x 1360 mm, 27" x 28" x 54" (L x D x H) *
Top Shelf Dimension*	640 mm x 450 mm, 25" x 18" (L x D) *	450 mm x 735 mm 18" x 29" (L x D) *
Other 4 Shelves*	510 mm x 450 mm, 20" x 18" (L x D) *	450 mm x 610 mm 18" x 24" (L x D) *
MS Drawer at Bottom*	510 mm x 450 mm x 195 mm, 20" x 18" x 8" (L x D x H)*	450 mm x 610 mm x 195 mm 18" x 24" x 8" (L x D x H)*

Power In

- Electrical shock protection – class 1
- Degree of shock protection – Type B
- Ingress protection – IPX4
- Two Inputs (One for one output with 32A, DP and other with 6 output connections with 20A) with safety fuse and switch
- Power On/Off switch available on each output with indicator

Standard Features

- Adjustable shelf, SS rails between top shelf and drawer
- Drawers are available with lock and key mechanism
- Suitable for keeping light source, camera unit, VCR, printer, stabilizer etc.
- Five MS shelves
- Aesthetically designed vertical uprights of aluminum extrusion
- Four 125 mm diameter single wheel high grade synthetic body castors, two with brake and without brake
- Stainless steel IV pole and gas cylinder cage
- Two stainless steel push handles

Add on Accessories



Camera Hanger



Keyboard Tray 9012, 9014



Pipe Spanner Hanger

- Tuber hanger in front top shelf

*±10mm Engineering variation



Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

We are interested in growing relationships



OTS-B:05:19.05.23