

OT Solutions



JANAK®
Because we care.



SURGISTAR
MULTI PURPOSE ELECTRO-MECHANICAL OPERATING TABLE



Provision for Ortho Attachment



User Friendly Handset



Railings with Locking System for Safety



Auxiliary Control on Column



Manual Override System Control



Twin Disk Castors with Central Locking Lever 2 Swivelling and 2 Fixed

Surgical procedures has become customizable with optimized workflow, improved efficacies yet cost-effective measures resulting in maximized profitability. Surgistar Electric Operating Table systems are unmatched in versatility, with detachable and interchangeable head and Leg sections for various surgical procedure, including radio translucent table top that offers outstanding C-arm imaging/ fluoroscopy. Surgistar tables can be adapted to the desired operating environment and are capable of communicating with C-Arm systems to provide the ultimate in imaging and surgical performance. With unmatched adaptive positioning options, Surgistar facilitates seamless patient positioning and improves ergonomics for the surgical team.

Surgistar provides state-of-the-art solutions for cutting-edge, configurable surgical table systems. Surgistar multifunctional OT Tables enable flexible OR function for both minimally invasive and open surgical procedures.



Trendelenburg / Reverse Trendelenburg



Lateral Tilt (Left / Right)



Height Up & Down



Lateral Thorax Position for Kidney and Thoracic Surgery

OT Solutions

Technical Specifications

Table Top Size (LxW)*	1995 x 500 mm
Minimum Height (Without Mattress)*	745 mm
Maximum Height (Without Mattress)*	995mm
Trendelenburg / Reverse Trendelenburg**	25° / 25°
Lateral Tilt**	20° / 20°
Head Section Adjustment**	+45° / -90°
Back Section Adjustment**	+75° / -25°
Leg Section Adjustment**	+20° / -90°
Kidney Bridge*	100 mm
Weight Bearing Capacity^	200 kgs

Standard Accessories

"L" Shaped Anaesthetic Screen	One pc
Padded Shoulder Support	One pair
Padded Lateral Support	One pair
Radio-Translucent Armboard with Pad	Two pcs
Goepel Knee Crutch with Straps	One pair
PU Foam Mattress	One set
Body Strap	One pc
Power supply	220V AC - 50Hz

Variant

Surgistar Core

4 Motorised Functions (Height, Trend/Rev Trend, Lateral Tilt and Back Section), Flex-Reflex Position, Kidney Bridge, Single Leg Section, Handset, Auxillary Control on Column, Battery Backup

Surgistar Pro

4 Motorised Functions (Height, Trend/Rev Trend, Lateral Tilt and Back Section), Kidney Bridge, Single Leg Section, Handset, Auxillary Control on Column, Battery Backup, '0' Position and Two Programmable Memory Function

Surgistar Advance

4 Motorised Functions (Height, Trend/Rev Trend, Lateral Tilt and Back Section), Kidney Bridge, Single Leg Section, Handset, Auxillary Control on Column, Flex-Reflex Position, Battery Backup, Hand Crank Operated Manual Override System


* ±10 mm Engineering variation. ** ± 2° Engineering variation. ^When operating table is configured with standard accessories

Additional Accessories for Operating Table **Surgistar**



For other accessories contact Midmark India sales representative

C-Arm systems to provide the ultimate in imaging and surgical performance



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

We are interested in growing relationships



Surgistar: 01:16.02.23