



1036- SS Rigid Wheelchair with Cushioned Seat and Back

- 01 Overall approximate dimensions: 1120 mm x 680 mm x 950 mm, 44" x 27" x 37" (L x W x H)*
- 02 SS tubular framework fitted with cushioned seat and back
- 03 Adjustable aluminum footrest
- 04 Two solid rubber tyred bicycle wheels with brake and self-propelling SS hoops
- 05 Two swivel castors in front of 200 mm diameter, 8"
- 06 SS tubular construction
- 07 Safety belt provided
- 08 Weight bearing capacity: 135 Kgs
- 09 Upholstery colour: Robust Brown

Add on Accessories

- 01 Moulded chart holder



*The product images shown are for illustrative purpose only and not be an exact representation of the product. Accessories shown in the product image may not be part of the standard configuration.

We are interested in growing relationship



1033A



1033A

1033A – MS Rigid Wheelchair with Cushioned Seat and Back

- 01 Overall approximate dimensions: 1120 mm x 680 mm x 950 mm, 44" x 27" x 37" (L x W x H)*
- 02 MS tubular framework fitted with cushioned seat and back
- 03 Adjustable aluminium footrest
- 04 Two solid rubber tyred bicycle wheels with brake and self-propelling SS hoops
- 05 Two swivel castors, 200 mm* diameter in front
- 06 Safety belt provided
- 07 Pre-treated and epoxy powder coated
- 08 Supplied in Knocked Down Condition (KDC)
- 09 Weight bearing capacity: 135 Kgs
- 10 Upholstery colour: Robust Brown

Add on Accessories

- 01 Moulded chart holder

*± 10 mm Engineering Variation and ** ± 5° Engineering Variation

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

1033A 1036:09:28.02.24