

Transfer Solutions

1040AR

Lightweight aluminum wheelchair
with 30° reclining backrest



1040AR



1040 AR is a lightweight aluminum wheelchair
with 30° reclining backrest and foldable to the front.

*The product images shown are for illustrative purpose only and not be an exact representation of the product.
Accessories shown in the product image may not be part of the standard configuration.

We are interested in growing relationship





01 30° backrest inclination



02 Padded nylon seat and backrest



03 Swing-away and detachable armrests



04 In height and depth adjustable arm pad



05 In and outwards swingable and detachable leg rest



06 Push Handles




07 Front wheel



08 Metal push bar



09 Integrated lowering of the rear axle



*± 10 mm Engineering Variation and ** ± 5° Engineering Variation

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
 Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

1040AR:00:05.04.23

