

Transfer Solutions

# 1040CR

Multi Comfort Wheelchair



1040CR

**1040CR** is a highly comfortable wheelchair with multiple adjustment possibilities making it adaptable to almost all body types. It has an adjustable backrest from 0° to 54° and an adjustable seat from 0° to 21°. Thanks to the pre-formed cushion on the seat and the back, its leg supports with comfortable pillows, and its fully adjustable headrest, the 1040 CR has been designed to offer the user maximum comfort.



01 Complete tilting for total relaxation

\*The product images shown are for illustrative purpose only and not be an exact representation of the product. Accessories shown in the product image may not be part of the standard configuration.

We are interested in growing relationship





02 Adjustable headrest in inclination (angle), height and depth with rotation and translation



03 Pushing handles adjustable in height



04 Side Support




05 Tilting footplates



06 Adjustable and removable leg rests



07 Ergonomically shaped memory foam cushion



08 Individual positioning through the adjustable backrest inclination



09 Adjustable tipping piece at the back



13 Easy to disassemble for transport or storage



10 Front wheel



11 16" Rear wheel



12 Additional accessories



Total Width  
680 mm



Seat Width  
460 mm




Dismantle Width  
600 mm



Inclination  
Backrest 0°-54°



Seat Height  
460-530 mm




Total Height  
990-1060 mm



Total Length  
1050-1090 mm



Seat Depth  
460-520 mm



Backrest Height  
550-585 mm




Users Weight  
135 Kg



Weight  
44.3 Kg



Obstacle Climbing  
60 mm



Turning Radius  
1520 mm



Rear Wheels  
16"

\*± 10 mm Engineering Variation and \*\* ± 5° Engineering Variation

**Midmark (India) Private Limited**

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070  
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

1040CR:00:02:05:23

