



1704

1704 – Stretcher on Trolley (Fully SS)

- Overall approximate dimensions: 2050 mm x 570 mm x 820 mm, 81" x 22" x 32" (L x W x H)*
- Mattress platform dimension: 1810 mm x 560 mm, 71" x 22" (L x W)*
- Mattress dimension: 1830 mm x 560 mm x 50 mm (thickness with two safety belts), 72" x 22" x 2"
- Safe working load of 135 Kgs and patient weight bearing capacity of 130 Kgs
- SS tubular framework mounted on four swivel castor of 200 mm, 8" diameter, two with brake and two without brake
- Removable stretcher top made of stainless steel sheet supported on SS tubular frame
- Two provisions for IV pole
- Supplied with safety belts



1703

1703- Stretcher on Trolley

- Overall approximate dimensions: 2030 mm x 570 mm x 820 mm, 80" x 22" x 32" (L x W x H)*
- Mattress platform dimension: 1810 mm x 560 mm, 71" x 22" (L x W)*
- Mattress dimension: 1830 mm x 560 mm x 50 mm (thickness with two safety belts), 72" x 22" x 2"
- Safe working load of 135 Kgs and patient weight bearing capacity of 130 Kgs
- Mild steel tubular framework mounted on four swivel castors, 200 mm (8") diameter, two with brake and two without brake
- Removable stretcher top made of aluminium sheet supported on tubular frame
- Pre-treated and powder coated
- Supplied in KDC
- Supplied with safety belt

Add on Accessories

- | | |
|--|--|
| <ul style="list-style-type: none"> • Rexine covered 50 mm, 2" thick single mattress 1830 mm x 560 mm, 72" x 22" (L x W) for 1704 & 1703 • SS IV Pole, fixed height for 1704 and 1703 • SS Oxygen cylinder cage for 1704 • MS Oxygen cylinder cage for 1703 | <ul style="list-style-type: none"> • MS storage tray for 1703 • SS storage tray for 1704 • Swing-away type railing with SS tubular frame for 1704 (Pair) • Swing-away type railing with MS tubular frame for 1703 (Pair) |
|--|--|





1702A – Stretcher on Trolley

- Overall approximate dimension: 1845 mm x 560 mm x 845 mm, 72" x 22" x 33" (L x W x H)*
- Mattress platform dimension : 1790 mm x 555 mm, 70" x 22" (L x W)*
- Mild steel tubular framework mounted on four swivel castors, 200 mm (8") diameter two with brake and two without brake
- Removable stretcher top in two section made of pre-laminated board supported on tubular frame
- Backrest on ratchet
- Supplied with safety belts
- Fitted with swing away type mild steel epoxy powder coated side rails 630 mm (635 mm with locking flat)
- Oxygen cylinder cage & MS storage tray
- Stainless steel IV poles with the provisions on diagonal sides
- Pre-treated and powder coated finish
- Supplied in KDC
- Safe working load of 135 Kgs and patient weight bearing capacity of 105 Kgs

Add on Accessories

- Rexine covered mattress of 50 mm (2") *thickness for 1702A



1712- Folding Stretcher Trolley

- Overall approximate dimension: 2100 mm x 640 mm, 83" x 25" (L x W)* (Without railing)
- Folded length approximately 1020 mm, 40" (L)*
- Fixed Height: 830 mm, 33"*
- MS tubular framework mounted on four swivel castors of 200* mm diameter, two with brake and two without brake
- Fixed foldable stretcher top made up of aluminium sheet powder coated bolted on tubular trolley framework
- Folding on ratchet at both ends
- Provision for IV poles
- Pre-treated and powder coated
- Supplied with safety belts
- Safe working load of 135 Kgs and patient Weight bearing capacity of 125 Kgs

Add on Accessories

- Rexine covered 50 mm (2") thick single section mattress 1850 mm x 560 mm, 73" x 22" (L x W)*
- IV pole for stretcher trolley
- Swing-down type railing with SS tubular frame (pair)

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.



* ±10 mm Engineering variation. ** ± 5° Engineering variation. ^When stretcher trolley is configured with standard accessories

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 www.midmark.in | india@midmark.com

TS-B:05:01.01.24

