



1707R at maximum height



1707R at minimum height

1707R- Emergency & Recovery Trolley

Technical Specifications

Overall Approximate Dimension with Moulded Handle*	2000 mm x 765 mm, 79" x 30" (L x W)
Overall Approximate Dimension*	1890 mm x 815 mm, 74" x 32" (L x W)
Stretcher Dimension*	1845 mm x 860 mm, 73" x 34" (L x W)
Mattress Platform Dimension *	1790 mm x 560 mm, 70" x 22" (L x W)
Backrest**	0° - 65°
Trendelenburg**	0° - 11°
Reverse Trendelenburg**	0° - 11°
Minimum Height*	625 mm, 25"
Maximum Height*	870 mm, 34"
Safe Working Load	165 Kgs
Patient Weight Bearing Capacity^	135 Kgs

Add on Accessories

- Removable and detachable monitor shelf
- Rexine covered mattress of 50 mm (2") thickness with two sections

Safety Standards

- Compliance for safety: International Standards IEC 60601-2-52

Standard Features

- Height adjustment by screw mechanism on chromium plated forged handle
- MS tubular framework
- Removable two section x-ray top with safety belts
- Sliding x-ray cassette holder
- Trendelenburg and reverse trendelenburg positions assisted by gas spring
- Raised backrest on ratchet at 25°, 35°, 45°, 55° & 65°
- MS storage tray
- 100 mm buffers at 4 corners
- High grade synthetic body castors (125 mm) two with brakes and two without brakes
- SS telescopic IV pole with two hooks along with provision at four positions
- Epoxy powder coated oxygen cylinder cage
- Collapsible side railing with MS tubular frame

CE *± 10 mm Engineering Variation ** ± 5° Engineering Variation ^ When stretcher trolley is configured with standard accessories

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

TS-A:00.19.12.23

We are interested in growing relationship



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.