



AstirONXT

AstirO NXT with SS fully collapsible safety side railing (pair)

Stretcher trolleys play an essential role in the mobility and transfer of patients during emergencies. AstirO NXT is engineered for quick and safe intrahospital patient transfer equipped with emergency and recovery care management. AstirO NXT comes with height adjustment on hydraulic columns. Stable and robust it is designed to intensive use, making it a perfect choice for critical care and high dependency areas.



AstirO NXT with big size polymer safety side railing (pair)

Technical Specifications

Overall Dimension*	2060 mm x 870 mm (L x W)
Stretcher Dimension*	2000 mm x 980 mm (L x W)
Mattress Platform Dimension*	1855 mm x 700 mm (L x W)
Backrest**	0° - 65°
Trendelenburg**	0° - 14°
Reverse Trendelenburg**	0° - 14°
Height*	630 mm to 950 mm (without mattress)
Safe Working Load	170 kgs
Patient Load Bearing Capacity^	150 kgs

Variants

- AstirO NXT with SS Railing
- AstirO NXT with Polymer Railing

Add on Accessories

- Rexin covered two section PU foam mattress, 50 mm thick
- SS heavy duty IV pole with 4 hooks

Standard Features

- Height adjustment by hydraulic pump with foot pedal on both sides
- Gas spring assisted manual backrest adjustment
- TR/RTR operation by hydraulic pump with foot pedal on both sides
- Interchangeable SS safety board and push handles on both head and leg end
- Two section mattress platform made of X-Ray permeable pre-laminated board, supported on strong tubular frame with sliding multisize X-Ray cassette holder
- SS heavy duty IV pole with 2 hooks
- 200 mm diameter single wheel castor with central locking mechanism, operational from both sides
- Safety corner buffers & push pull handles with rubber grippers (Available in red & blue colours)
- Oxygen cylinder cage
- Safety belts, pair

Safety Standard

Compliance for safety:
International standards IEC 60601-2-52



Maximum Height



Minimum Height



Trendelenburg



Reverse Trendelenburg



Oxygen Cylinder Cage



X-Ray Permeable Platform



Large Diameter Castors
with Multistage Lock



Safety Head Board



Safety Side Railings



Safety corner buffers & push pull handles with rubber grippers
(Available in red & blue colours)



* ±10 mm Engineering variation. ** ± 5° Engineering variation. ^When stretcher trolley is configured with standard accessories

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

AstirO NXT: 06:08.12.23

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.