

Transfer Solutions



JANAK[®]
Because we care.



1761

The safe way to transfer patients

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt.Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalog in any form is liable for legal action.

Transfer Solutions



Midmark Transfer Systems

make day-to-day work with
patients easier & efficient

The safe way to transfer patients

In a hospital, the working conditions should be effortless and comfortable. Caregivers should avoid lifting the patients while transferring.

Midmark Rollbord™ assists the caregiver to transfer the patients without lifting, ensuring comfort and safety of the patients. These products are designed to facilitate the transfer without using a lot of effort, thus making the work easier and efficient for the caregiver's and reduces stress on their back, thus ensuring efficiency within critical care units.

Midmark Transfer Systems — a great help to healthcare professionals

- Prevents back injuries to caregiver's
- Skin friendly (no latex)
- Light weight
- Easy to clean and disinfect— hygiene approved



Turn the patient onto the side, slide the board under the body and return the patient to the original position



Gently push the patient at the shoulder and hip, without applying force, and transfer.



Turn the patient onto the side. Remove the Rollbord™ starting at the feet.

Transfer Solutions

The hightec Rollbord™

To suit every clinical application

Standard Foldable



1761

Suitable for radiology, operating theatre, intensive care units and all special wards.

Size: 177x50cm/89x 50cm, 2.7 kg 70x20 inch/35 x20 inch

Safe Working load: 200 Kgs

- Anti-bacteriologic
- Anti-fungic

Environmentally friendly material*

- MIDMARK transfer aids are made of ecologically safe materials.
- Inner Rollbord™ core is made of plastic, coated with glass fibre fabric. Roll material is made of tear resistant & antibacterial nylon fabric, assisting in safe and smooth transfer.
- Teflonised sliding fabric around a polyethylene foam core, coated polyester cloth, PE security endings, which makes sliding easy.

Easy Cleaning

- The Rollbord™ can be cleaned best with wet wipes or any hospital disinfectant.

Transfer Solutions

Features which suit every hospital application —

- High performance
- Radioluscent
- Non-toxic materials
- Anti-bacterial surface
- Safe and comfortable for patients
- Prevents back injuries to caregiver's
- Skin-friendly (no latex)
- Easy to clean and disinfect — hygiene approved
- Disposable covers available
- Tear resistant nylon cloth
- Bridges gaps up to 20 cm
- Lightweight & easily suspendable

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt.Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalog in any form is liable for legal action.

Midmark (India) Private Limited

Head Office: Art Guild House, A-Wing, Unit No. 15 & 16, Phoenix Market City, LBS Road, Kurla (W), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com
Check the nearest **Midmark India Experience Centre** at www.midmark.in/miec

We are interested in growing relationships



TS-F;01;08.08.19