

Transfer Solutions

MX4000B

Multi Treatment Chair



MX4000B

Transforming care experiences with MX4000B

In the realm of healthcare, comfort plays a pivotal role in providing care. Designed to address your and your patient's needs in mind, MX4000B empowers you to deliver impeccable care with precision. Your path to unrivaled patient care excellence.

Space Optimization:

In outpatient and ambulatory care settings where space is limited, MX4000B stands out with its compact design, maximizing the efficiency of your facility's space.

Elevated Aesthetics:

MX4000B's visually appealing design adds a touch of warmth to healthcare environments, uplifting the spirits of both patients and caregivers.

Patient Ease and Wellbeing:

Navigating treatments such as chemotherapy or dialysis can be tough. The environment is key. Swapping a hospital bed for a less intimidating treatment chair can lessen emotional stress, allowing patients to concentrate on recovery. This small but meaningful change can greatly improve a patient's overall well-being.

Seamless Mobility and Accessibility:

With its user-friendly design, MX4000B can be effortlessly repositioned or transported across different sections of the healthcare facility.



*The product images shown are for illustrative purpose only and not be an exact representation of the product. Accessories shown in the product image may not be part of the standard configuration.

We are interested in growing relationship




Patient Positioning. Comfort. Accessibility.



Supine Position



Trendelenburg



Extendable Foot Support



Reclining Position



Chair



Handset for Easy Positioning

Technical Specifications

| | |
|--------------------|---|
| Overall Dimension* | 1860 mm x 850 mm, 73.22" x 33.46" (LxW) |
| Seat Width* | 600 mm, 23.62" |
| Extended Length* | 1920mm, 75.59" |
| Minimum Height* | 620mm, 24" |
| Maximum Height* | 740mm, 29" |

Electrical Parameters

| | |
|-------------------------------|--|
| Power In (Mains Power Supply) | Nominal 240 VAC |
| Operating Range with Battery | 90 VAC to 300 VAC 45/55 Hz, Max 2A |
| Protection Voltage | 100 VAC to 240 VAC |
| Electrical Shock Protection | Class I |
| Degree of Shock Protection | Type B |
| Liquid Ingress Protection | IPX6 |
| Rechargeable Batteries | 2 x 12 volt Sealed Lead/Acid Gel |
| Duty Cycle | 10 % (Two minutes for every eighteen minutes) |

Standard Accessories



Paper Roll Holder

Add-On Accessories



Food Table / Reading Table



IV Pole



Safety Belts

Upholstery Colour Options



Raisin Grape + Navy



Healing Waters + Robust Brown

Features (Inbuilt)

| | |
|---|--|
| Adjustable Armrests with Back & Forth and Up-Down Movements | Available |
| Patient Hand Control | Available |
| Chair Position | Available |
| Emergency Head Down | Available |
| Castors* | 125 mm high grade synthetic, single wheel castors, four with brake |
| Backrest** | 0° - 72° |
| Trendelenburg** | 0° - 15° |
| Battery Back-Up | Available |
| Safe Working Load | 150 kgs |
| Patient Load Bearing Capacity^ | 135 kgs |

Safety Standards

| | |
|---|---------------|
| General safety for medical electrical equipment | IEC 60601-1-1 |
| EMI/EMC | IEC 60601-1-2 |

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.

* ±10 mm Engineering variation. ** ± 5° Engineering variation. ^When bed is configured with standard accessories. #Factory fitted.

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

MX4000B:03:17.11.23

