

Wardcare Solutions



JANAK[®]
Because we care.



1043A

1043A – Stainless Steel Saline Stand with Two Hooks (Tubular Base)

- Strong stainless steel tubular construction mounted on five pronged rectangular tubular base
- Fitted with five swivel, rust proof, synthetic body castors 50 mm (2") wheel diameter without brake
- Stainless steel IV Pole with two hooks
- Height adjustable from approximate 1395 mm (min)* to 2290 mm (max)*
- Load capacity: 10 Kgs



1044

1044- Mild Steel Saline Stand with Two Hooks (Tubular Base – Powder Coated)

- Strong mild steel tubular construction mounted on five pronged rectangular tubular base
- Fitted with five swivel, rust proof, synthetic body castors 50 mm (2") wheel diameter without brake
- Stainless steel IV pole with two hooks
- Height adjustable from approximate 1395 mm (min)* to 2290 mm (max)*
- Pre-treated and powder coated finish
- Load capacity: 20 Kgs

*± 10 mm Engineering Variation



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.



Wardcare Solutions



1043B

1043B - Stainless Steel Saline Stand with Four Hooks (Tubular Base)

- Strong stainless steel tubular construction mounted on five pronged rectangular tubular base
- Fitted with five swivel, rust proof, synthetic body castors 50 mm (2") wheel diameter without brake
- Stainless steel IV pole with four hooks
- Height adjustable from approximate 1395 mm (min)* to 2290 mm (max)*
- Load capacity: 20 Kgs



1044A

1044A - Mild Steel Saline Stand with Four Hooks (Tubular Base – Powder Coated)

- Strong mild steel tubular construction mounted on five pronged rectangular tubular base
- Fitted with five swivel, rust proof, synthetic body castors 50 mm (2") wheel diameter without brake
- Stainless steel IV pole with four hooks
- Height adjustable from approximate 1395 mm (min)* to 2290 mm (max)*
- Pre-treated and powder coated finish
- Load capacity: 20 Kgs

*± 10 mm Engineering Variation



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

We are interested in growing relationships



WCS-E:06:19.05.23