Wardcare Solutions





AX1016A Deluxe Adjustable Bedside Table

- Overall size: 810 mm x 355 mm, 32 " x 14" (L x W)*
- Two section laminated top with railing on fixed section.
 Raised edge on one side on tiltable section
- Aluminium extruded telescopic section fitted with gas spring for height adjustable from approximately 815 mm to 1125 mm, 32" to 44"
- Mild steel rectangular tubular base frame mounted on four swivel synthetic body castors of 50 mm
- · Pre-treated and powder coated finish
- Supplied in KDC



AX1016P Deluxe Adjustable Bedside Table

- Overall size: 830 mm x 406 mm, 33 " x 16" (L x W)*
- Single section non adjustable polymer top.
- Aluminium extruded telescopic section fitted with gas spring for height adjustable from approximately 815 mm to 1125 mm, 32" to 44"
- Mild steel rectangular tubular base frame mounted on four swivel synthetic body castors of 50 mm
- · Pre-treated and powder coated finish
- Supplied in KDC

*± 10 mm Engineering Variation









 ϵ

Wardcare Solutions



AX1012A Wardcare Adjustable Bedside Table

- Overall dimension: 960 mm x 420 mm, 38" x 17" (L x W)* Top size: 810 mm x 355 mm, 32" x 14" (L x W)*
- Two section laminated top with railing on fixed section raised edge on one side tiltable section
- Mild steel tubular telescopic stem, with geared SS handle for height adjustment from approximately 715 mm to 1080 mm, 28" to 43"*
- · Mild steel rectangular tubular base frame mounted on four twin wheel non rusting castors of 50 mm
- Pre-treated and powder coated
- Supplied in KDC



AX1012P Wardcare Adjustable Bedside Table

- Overall dimension: 910 mm x 406 mm, 36" x 16" (L x W)* Top size: 810 mm x 405 mm, 32" x 16" (L x W)*
- Single section non adjustable polymer top.
- · Mild steel tubular telescopic stem, with geared SS handle for height adjustment from approximately 715 mm to 1080 mm, 28" to 43"*
- · Mild steel rectangular tubular base frame mounted on four twin wheel non rusting castors of 50 mm
- Pre-treated and powder coated
- Supplied in KDC



AX1011 Wardcare Overbed Table

- Overall approximate dimension: 1120 mm x 450 mm x 1070 mm, 44" x 18" x 42" (L x W x H)*
- Laminated top of approximate 760 mm x 450 mm, 30" x 18" (L x W)*
- Mild steel tubular base frame mounted on four swivel synthetic body castors of 50 mm
- · Pre-treated and powder coated finish
- Supplied in KDC

*± 10 mm Engineering Variation

CE





