

# Wardcare Solutions



**JANAK®**  
Because we care.



**AX1016A**

## AX1016A Deluxe Adjustable Bedside Table

- Overall size: 810 mm x 355 mm, 32" x 14" (L x W)\*
- Two section laminated top with railing on fixed section. Raised edge on one side on tiltable section
- Aluminium extruded telescopic section fitted with gas spring for height adjustable from approximately 815 mm to 1125 mm, 32" to 44"
- Mild steel rectangular tubular base frame mounted on four swivel castors
- Pre-treated and powder coated finish
- Supplied in KDC

Tabletop Colour options :  Sky Blue  Peppermint  Wheat Oak  New European Ash



**AX1016C**

## AX1016C Elite Bed Side Table

- Overall size: 750 mm x 360 mm, 30" x 14" (LxW)\*
- Single section nonadjustable laminated top with SS railing on three sides.
- Telescopic section fitted with gas spring for height adjustable from approximately 785mm to 1115mm, 31" to 45"
- Mild steel round base frame mounted on four swivel castors
- Pre-treated and powder coated finish
- Supplied in KDC

Tabletop Colour options :  Sky Blue  Peppermint  Wheat Oak  New European Ash  Metallic White



**AX1016P**

## AX1016P Deluxe Adjustable Bedside Table

- Overall size: 910 mm x 405 mm, 36" x 16" (L x W)\*
- Single section non adjustable polymer top.
- Aluminium extruded telescopic section fitted with gas spring for height adjustable from approximately 835 mm to 1150 mm, 33" to 45"
- Mild steel rectangular tubular base frame mounted on four swivel castors
- Pre-treated and powder coated finish
- Supplied in KDC

\*± 10 mm Engineering Variation



# Wardcare Solutions



AX1012A

## AX1012A Wardcare Adjustable Bedside Table

- Overall dimension: 960 mm x 420 mm, 38" x 17" (L x W)\*  
Top size: 810 mm x 355 mm, 32" x 14" (L x W)\*
- Two section laminated top with railing on fixed section raised edge on one side tiltable section
- Mild steel tubular telescopic stem, with geared SS handle for height adjustment from approximately 715 mm to 1080 mm, 28" to 43"\*
- Mild steel rectangular tubular base frame mounted on four twin wheel castors
- Pre-treated and powder coated
- Supplied in KDC

Tabletop Colour options : ● Sky Blue ● Peppermint ● Wheat Oak ● New European Ash



AX1012P

## AX1012P Wardcare Adjustable Bedside Table

- Overall dimension: 910 mm x 406 mm, 36" x 16" (L x W)\*  
Top size: 810 mm x 405 mm, 32" x 16" (L x W)\*
- Single section non adjustable polymer top.
- Mild steel tubular telescopic stem, with geared SS handle for height adjustment from approximately 730 mm to 1090 mm, 29" to 43"\*
- Mild steel rectangular tubular base frame mounted on four twin wheel castors
- Pre-treated and powder coated
- Supplied in KDC



AX1011

## AX1011 Wardcare Overbed Table

- Overall approximate dimension: 1120 mm x 450 mm x 1070 mm,  
44" x 18" x 42" (L x W x H)\*
- Laminated top of approximate 760 mm x 450 mm, 30" x 18" (L x W)\*
- Mild steel tubular base frame mounted on four swivel castors
- Pre-treated and powder coated finish
- Supplied in KDC

Tabletop Colour options : ● Sky Blue ● Peppermint ● Wheat Oak ● New European Ash

\*± 10 mm Engineering Variation



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.

**Midmark (India) Private Limited**

**Head Office:** Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070  
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com



WCS:09:24.08.23