

# Wardcare Solutions



**JANAK**<sup>®</sup>  
Because we care.



## BX1000 Simple Bed

The BX1000 simple bed is designed for remote primary healthcare centre, where beds are expected to last for years without any maintenance. This basic bed is specially designed to be strong, functional, and effective using the best quality raw materials and processes. Metal treatment and powder coating ensures the excellent finish.

## Technical Specifications

	BX1000	BX1000A
Overall Approximate Dimension*	2100 mm x 910 mm x 600 mm, 83" x 36" x 24" (L x W x H)	1940 mm x 760 mm x 600 mm, 76" x 30" x 24" (L x W x H)
Bed Frame Dimension*	1980 mm x 910 mm, 78" x 36" (L x W)	1830 mm x 760 mm, 72" x 30" (L x W)
Mattress Dimension*	1980 mm x 910 mm, 78" x 36" (L x W)	1830 mm x 760 mm, 72" x 30" (L x W)
Fixed Height *	600 mm, 24"	600 mm, 24"

## Standard Features

- Pre-treated and powder coated finish
- Sockets provided for mosquito net poles
- Safe Working Load is 150 kg
- Patient Working Load is 140 kg

## Add-On Accessories

- Mattress arrester (set of two)
- Single mattress with PU foam covered with PVC rexine
- Telescopic IV pole with suitable bracket
- Mosquito net poles (set of four)
- MS full length safety side railings epoxy powder coated for BX1000
- Crib with attachment for BX1000 series
- 3" PU foam with PVC rexine mattress for crib with attachments
- Provision for urine bag holder on either side of bed

\*± 10 mm Engineering Variation and \*\* ± 5° Engineering Variation



# Wardcare Solutions



## BX1001 Attendant Bed

The attendant bed is designed for remote primary healthcare centre, where beds are expected to last for years without any maintenance. This bed is specially designed to be strong, functional, and effective using the best quality raw materials and processes

### Technical Specifications

Code No.	Overall & Bed Frame Approximate Dimension in mm			
	L	W	H	(L x W)
BX1001	1980	910	400	78" x 36"
BX1001A	1830	762	400	72" x 30"
BX1001B	1830	838	400	72" x 33"
BX1001C	1905	762	400	75" x 30"

Mattress Dimension*	
BX1000 and BX1001	1980 mm x 910 mm, 78" x 36" (L x W)*
BX1000A and BX1001A	1830 mm x 762 mm, 72" x 30" (L x W)*
BX1001B	1830 mm x 838 mm, 72" x 33" (L x W)*
BX1000A	1830 mm x 760 mm, 72" x 30" (L x W)*
BX1001C	1905 mm x 762 mm, 75" x 30" (L x W)*

### Standard Features

- Safe working load -150 Kgs
- Patient working load -140 kgs
- Single section mattress covered with PVC rexine: 75 mm, 3", 32 density

### Add-On Accessories

- Mattress arrestor (set of two)
- Provision for urine bag holder on either side of bed


\*± 10 mm Engineering Variation



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.

## Midmark (India) Private Limited

**Head Office:** Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070  
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com



We are interested in growing relationships

BX1000-BX1001-07:19.05.23