



BX1003 Attendant Bed Cum Chair

Technical Specifications

Overall Approximate Dimension	
Dimension For End to End*	1850 mm x 625 mm, 73" x 25" (L x W)
Platform Dimension* (Without Mattress)	1850 mm x 635 mm, 73" x 25" (L x W)
Fixed Height*	515 mm, 20" H (Without Upholstery)
Upholstery Thickness*	75 mm, 3" H
Backrest Angle**	0° - 70°
Safe Working Load	150 Kgs
Patient Load Bearing Capacity	150 Kgs

Standard Features

- Manual chair position
- Manual adjustment for back section and leg section
- Two armrest pivoted with backrest on each side
- Tray at the bottom
- All MS parts are pre-treated & epoxy powder coated
- Upholstery colour - Navy Blue



0155B Crib with Attachment (on Clamp Provision)

- Overall crib size: 785 mm x 380 mm x 300 mm (30" x 15" x 12") (L x W x H)*
- Swings on two tubular rods fixed towards the leg side board of the bed
- Pre-treated and powder coated
- 0155 - crib with attachment for BX1000 series (specifications as per 0155B)
- 0122B - 3" PU foam with PVC rexine mattress



1051 Bed Side Screen (Three Fold)

- Overall size: 2440 mm x 1690 mm (96" x 67") (W x H)*
- Mild steel tubular construction in three sections
- Mounted on six swivel twin wheel non rusting castors, 50 mm* (2") diameter
- Middle frame 1220 mm* (48") wide
- Side spans of 610 mm* (24") each frame with four hooks
- Supplied with hooks and springs with curtains
- Pre-treated and powder coated
- Supplied in KDC



1069 - Wardcare Cylinder Trolley (Push Type)

- Suitable for 1320 litre gas cylinder
- Mild steel tubular framework mounted on 125 mm* (5") diameter castors
- Trolley provided with locking provision and stainless steel base
- Pre-treated and powder coated

CE *± 10 mm Engineering Variation

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 www.midmark.in | india@midmark.com

WCS-B:10:29.02.24



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.