

# Wardcare Solutions



**JANAK**<sup>®</sup>

Because we care.



## BX1100 Isolation Bed

The BX1100 heavy duty isolation bed with commercial ply top and MS integrated head & foot panel is specially designed to be strong, functional and effective using the best quality raw materials and processes.

### Technical Specifications

Overall Approximate Dimension	1850 mm x 940 mm x 415 mm, 73" x 37" x 16" (L x W x H)*
Bed Frame Dimension	1850 mm x 940 mm, 73" x 37" (L x W)*
Mattress Dimension	1780 mm x 895 mm, 70" x 35" (L x W)*

### Add-On Accessories

Single section mattress with 3" thick PU foam

### Standard Feature

• Powder-coated finish

\*± 10 mm Engineering Variation

#### Midmark (India) Private Limited

**Head Office:** Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070  
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | [www.midmark.in](http://www.midmark.in) | [india@midmark.com](mailto:india@midmark.com)

We are interested in growing relationships



BX1100:00:20.07.21