

# BX2000 Plus

## Semi Fowler Bed



Basic back adjustable bed with conventional head and foot boards and side railing



Basic back adjustable bed with stainless steel laminated head and foot boards and side railing



Basic back adjustable bed with polymer moulded head and foot boards, and side railing

## BX2000 Plus Semi Fowler Bed

The BX2000 plus semi fowler bed has the facility to raise the back sections, which is operated by the crank mechanism provided at the foot end of the bed. The raising section allows the nursing staff to change positions for various activities and clinical positions with relative ease.



\*The product images shown are for illustrative purpose only and not be an exact representation of the product. Accessories shown in the product image may not be part of the standard configuration.

We are interested in growing relationship



## Technical Specifications

Bed Frame Dimension*	1980 mm x 920 mm, 78" x 36" (L x W)
Mattress Platform*	1880 mm x 810 mm, 74" x 32" (L x W)
Buffer to Handle Dimension*	2180 mm 86"
Buffer to Buffer* (when castors provided)	2040 mm x 970 mm, 80" x 38" (L x W)
Fixed Height	580 mm, 23"
Backrest Angle**	0° - 75°
Safe Working Load	150 kg
Patient Weight Bearing Capacity^	110 kg

## Add-On Accessories

### Technical Specifications

- Collapsible type epoxy powder coated vertical safety side railing (pair)
- H' type legs 910 mm (36") width & 600 mm (24") height with rubber shoes
- H' type legs 910 mm (36") width & 450 mm (18") height with rubber shoes
- H' type legs 910 mm (36") width & 600 mm (24") height with 125 mm high grade synthetic castors, two with brake and two without brake
- H' type legs 910 mm (36") width & 450 mm (18") height with 125 mm high grade synthetic castors, two with brake and two without brake

### Head and Foot Boards

- Polymer moulded head and foot boards
- SS laminated head and foot boards
- Conventional head and foot boards

## Standard Features

- Raised backrest
- Two section metallic perforated top
- Pre-treated and epoxy powder coated MS parts
- Provision for IV pole
- Removable handle
- Easily removable head and foot boards
- Supplied in SKD condition

## Safety Standards

General Safety for Medical Beds IEC 60601-2-52

- Two section mattress with 4" thick PU foam of 40 density covered with PVC
- Two section mattress with 3" thick PU foam of 32" density covered with PVC
- Telescopic IV pole with 2 hooks
- Telescopic IV pole with 4 hooks
- Orthopaedic traction pulley for trauma cases
- Fixed polymer moulded handle
- Removable handle
- Moulded chart holder
- Mattress arrestor (Set of 2)
- 4" safety rubber buffer at four corners
- Urine bag holder

\*± 10 mm Engineering Variation and \*\* ± 5° Engineering Variation ^ When Bed is configured with standard accessories  
Note: Option of fixed polymer moulded handle available in lieu of removable handle.

### Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070  
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 www.midmark.in | india@midmark.com

BX2000 Plus:09:24.12.23



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.