

# BX2000D



BX2000D Basic Bed

## BX2000D Basic Bed

The BX2000D Basic Bed has the facility to raise the back section which is operated by a crank mechanism provided at the foot end of the Bed. The raising section allow the nursing staff to change positions for various activities and clinical positions with relative ease.



\*The product images shown are for illustrative purpose only and not be an exact representation of the product. Accessories shown in the product image may not be part of the standard configuration.

We are interested in growing relationship



## Technical Specifications

Bed Frame Dimensions*	2005 mm x 920 mm, 79"x 36" (L x W)
Mattress Platform*	1880 mm x 810 mm, 74"x 32" (L x W)
Buffer to Buffer (when castors provided)*	2050 mm x 970 mm, 84"x 39" (L x W)
Buffer to Handle (L)*	2180 mm, 86"
Backrest Angle** (Motorized)	0° - 70°
Safe Working Load	150 kgs
Patient Weight Bearing Capacity <sup>^</sup>	110 kgs

## Add-On Accessories

### Railings and Legs

<ul style="list-style-type: none"> <li>Collapsible type epoxy powder coated railing (pair)</li> </ul>
<ul style="list-style-type: none"> <li>Collapsible type epoxy powder coated MS safety side railing with SS tubes (pair)</li> </ul>
<ul style="list-style-type: none"> <li>Swing down type epoxy powder coated MS safety side railing (pair) (factory fitted)</li> </ul>
<ul style="list-style-type: none"> <li>Collapsible type epoxy powder coated vertical safety side railing (pair)</li> </ul>
<ul style="list-style-type: none"> <li>Collapsible type full-length type SS railing with epoxy powder coated frame (pair)</li> </ul>
<ul style="list-style-type: none"> <li>Single set big size polymer railing (factory fitted)</li> </ul>
<ul style="list-style-type: none"> <li>Split type polymer moulded safety side railing (pair) (factory fitted)</li> </ul>
<ul style="list-style-type: none"> <li>H' type legs 910 mm (36") width &amp; 600 mm (24") height with rubber shoes</li> </ul>
<ul style="list-style-type: none"> <li>H' type legs 910 mm (36") width &amp; 450 mm (18") height with rubber shoes</li> </ul>
<ul style="list-style-type: none"> <li>H' type legs 910 mm (36") width &amp; 600 mm (24") height with 125 mm high grade synthetic castors, two with brake and two without brake</li> </ul>
<ul style="list-style-type: none"> <li>H' type legs 910 mm (36") width &amp; 450 mm (18") height with 125 mm high grade synthetic castors, two with brake and two without brake</li> </ul>

## Standard Features

<ul style="list-style-type: none"> <li>Motorized backrest</li> </ul>
<ul style="list-style-type: none"> <li>Two section metallic perforated top</li> </ul>
<ul style="list-style-type: none"> <li>Pre-treated and epoxy powder coated MS parts</li> </ul>
<ul style="list-style-type: none"> <li>Provision for IV pole</li> </ul>
<ul style="list-style-type: none"> <li>Easily removable head and foot boards</li> </ul>
<ul style="list-style-type: none"> <li>Supplied in SKD condition</li> </ul>

## Safety Standards

General Safety for Medical Beds	IEC 60601-2-52
---------------------------------	----------------

<ul style="list-style-type: none"> <li>Two section mattress with 4" thick PU foam of 40 density covered with PVC</li> </ul>
<ul style="list-style-type: none"> <li>Two section mattress with 3" thick PU foam of 32 density covered with PVC</li> </ul>
<ul style="list-style-type: none"> <li>Telescopic IV pole with 2 hooks</li> </ul>
<ul style="list-style-type: none"> <li>Telescopic IV pole with 4 hooks</li> </ul>
<ul style="list-style-type: none"> <li>Orthopaedic traction pulley for trauma cases</li> </ul>
<ul style="list-style-type: none"> <li>Fixed polymer moulded handle</li> </ul>
<ul style="list-style-type: none"> <li>Moulded chart holder</li> </ul>
<ul style="list-style-type: none"> <li>Mattress arrestor (set of 2)</li> </ul>
<ul style="list-style-type: none"> <li>Mosquito net poles</li> </ul>
<ul style="list-style-type: none"> <li>4" safety rubber buffer at four corners</li> </ul>
<ul style="list-style-type: none"> <li>Urine bag holder</li> </ul>

## Head and Foot Boards

<ul style="list-style-type: none"> <li>Polymer moulded head and foot boards</li> </ul>
<ul style="list-style-type: none"> <li>SS laminated head and foot boards</li> </ul>
<ul style="list-style-type: none"> <li>Conventional head and foot boards</li> </ul>

\*± 10 mm Engineering Variation and \*\* ± 5° Engineering Variation ^ When Bed is configured with standard accessories

Note: Option of fixed polymer moulded handle available in lieu of removable handle.

### Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070  
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 | Fax: +91 22 4915 3100 | www.midmark.in | india@midmark.com

BX2000D:02:04.12.23



Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.