


ICU Solutions

EX4411



EX4411 CS

The EX4411 hospital bed is designed and built integrating state of the art technologies to enhance the patient's comfort and improve caregiver's efficiency. Its outstanding performance is complemented by its ergonomic features and best-in-class quality, which complies with the highest global safety standards.



*The product images shown are for illustrative purpose only and not be an exact representation of the product. Accessories shown in the product image may not be part of the standard configuration.

We are interested in growing relationship



Technical Specifications

Overall Dimension*	2150 mm x 1025 mm (L x W)
Platform Frame*	2070 mm x 955 mm (L x W)
Mattress Platform*	1890 mm x 865 mm (L x W)
Mattress Platform Height* (Without Mattress)	395 mm to 695 mm
Backrest**	0° - 70°
Kneerest**	0° - 45°
Trendelenburg**	0° - 11°
Reverse Trendelenburg**	0° - 11°
Safe Working Load**	200 Kgs
Patient Load Bearing Capacity**	175 Kgs

Electrical Parameters

Power In (Mains Power Supply)	100 - 240 V ~, 50/60 Hz., Max 5A
Electrical Shock Protection	Class I
Degree of Shock Protection	Type B
Liquid Ingress Protection	IPX6
Duty Cycle	10 % (Two minutes for every eighteen minutes)

Standard Features

- One touch key for cardiac chair, auto CPR
- Motorized backrest, kneerest, height adjustment, trendelenburg & reverse trendelenburg
- Wired Handset control for patient and attendant.
- Dual side manually operated CPR
- Angle indicator for backrest, trendelenburg & reverse trendelenburg
- Full length polymer moulded split type safety side rails
- Polymer moulded head & foot panels, easily removable without using any tools
- Four corner safety rubber buffers
- Four section metallic perforated top
- 125 mm diameter single wheel high grade synthetic body castors with diagonal locking or central locking mechanism
- Inbuilt SMPS
- 100 mm thick four section mattress

Safety Standards

General safety for medical beds	IEC 60601-2-52
General safety for medical electrical equipment	IEC 60601-1-1
EMI/EMC	IEC 60601-1-2

Accessories and color shown in the picture may not be part of standard configuration. Technical alterations reserved in view of our continuous improvements. Midmark India Pvt. Ltd. is the licensed user of the trademark "Midmark India" and the logo. Reproduction of this catalogue in any form is liable for legal action.


• All variants are available with X-Ray permeable backrest (XR) feature which are factory fitted.

* ±10 mm Engineering variation. ** ± 5° Engineering variation. ^When bed is configured with standard accessories. #Factory fitted.


Add-On Accessories



Heavy Duty and Telescopic IV Pole



Vertical O2 Cylinder Cage



Orthopaedic Traction Pulley



Heavy Duty SS Lifting Pole



Moulded Chart Holder



Urine Bag Holder



Battery Backup#

Midmark (India) Private Limited

Head Office: Art Guild House, FOF, B-16, "B" Wing, 4th floor, Phoenix Market City, LBS Road, Kurla (West), Mumbai - 400 070
Toll Free: 1800 22 8020 | Tel: +91 22 4915 3000 www.midmark.in | india@midmark.com

EX4411: 05:25.09.23

